icare

tonometer VET

USER'S AND MAINTENANCE MANUAL

ENGLISH

FRANÇAIS

DEUTSCH

ESPAÑOL

TONOMETER

TonoVet tonometer TV01

INSTRUCTION MANUAL v.1.2 11/08

EN

CE 0044

Copyright © 2008 Icare Finland Oy
Made in Finland

Icare Finland Oy

Hevosenkentä 3, FI-02600 Espoo, Finland
Tel. +358 9 8775 1150, Fax +358 9 728 6670
Internet: www.icarefinland.com, e-mail: info_at_icarefinland.com

Contents

| | |
|---|----|
| Introduction..... | 3 |
| Safety instructions..... | 4 |
| Parts of the tonometer..... | 5 |
| Turning the device ON..... | 5 |
| Loading the probe..... | 6 |
| Turning the tonometer OFF..... | 6 |
| Measurement..... | 7 |
| Display after the measurements..... | 7 |
| Accessing old measurement values..... | 8 |
| Changing the calibration table for different species..... | 8 |
| Replacing the probe base..... | 9 |
| Replacing the batteries..... | 9 |
| Error messages..... | 10 |
| Service..... | 10 |
| Technical information..... | 11 |
| Spare parts and supplies..... | 11 |
| Diagram of tonometer functions..... | 11 |

Introduction

The TonoVet tonometer is used in the diagnosis, follow up, and screening of glaucoma in veterinary medicine (dogs, cats, horses etc.)

The TonoVet tonometer is based on the new patented induction based rebound method, which allows intraocular pressure (IOP) to be measured accurately, rapidly and without a local anaesthetic.

There is no risk of microbiological contamination, as single-use probes are used in the measurement.

Intraocular pressure changes due to the effects of pulse, breathing, eye movements, and body position. Because the measurement is made handheld in fractions of a second, several measurements are needed to obtain an accurate reading. The software is pre-programmed for six measurements, but the result can be seen from the first successful measurement.

There are calibration tables for different species (dog/cat, horse).


For veterinary use only.

Safety instructions

Read this manual carefully, because it contains important information on using and servicing the tonometer.

- Keep this manual for future use.
- When you have opened the package, check for any external damage or faults. Check especially that the casing is not damaged. If you suspect that there is something wrong with the tonometer, contact the manufacturer or the distributor.
- Use the tonometer only for measuring intraocular pressure. Any other use is incorrect. The manufacturer cannot be held responsible for any damage arising from improper use, or for the consequences thereof.
- Never open the casing of the tonometer, except for the battery compartment or to change the probe base.
- The manual contains instructions for replacing batteries and changing the probe base.
- Never use the tonometer in wet places.
- The probe base, battery compartment cover, collar and probes are so small that a child or animal could swallow them. Keep the tonometer out of the reach of children or animals.
- Do not use the device near the substances that are flammable - including flammable anaesthetic agents.
- Check that a clean new disposable probe is used before each measurement.
- Be sure that the probe contains the small plastic round tip in front.
- To avoid transmitting certain microbiological agents (e.g. bacteria), the device should be cleaned regularly with a disinfectant, e.g. an alcohol solution.
- The tonometer must not come into contact with the eye, except for the probes for a fraction of a second during the measurement, therefore do not bring the tonometer into contact with the eye or push it to the eye (tip of the probe should be 4–8 mm or 1/6 - 1/3 inch, from the eye).
- The tonometer conforms to EMC requirements (IEC 60101-1-2: 2001), but interference may occur in it if used near to (<1m) a device causing high-intensity electro magnetic emissions (like a cellular phone). Though the tonometer's own electromagnetic emissions are well below the levels in the relevant standards, they may cause interference in another device, e.g. a sensitive sensor, nearby.
- If the device is not to be used for a long time, removing the batteries is recommended, as the AA batteries may leak. Removing the batteries will not affect the subsequent functioning of the tonometer.
- Be sure to dispose of the single-use probes properly (e.g., in a container for disposable needles), because they may contain micro-organisms from the measured patient.
- Batteries, packaging materials, and probe bases must be disposed of according to local regulations.

Parts of the tonometer




Turning the device ON

Place the wrist strap into the wrist strap attachment. Place the wrist strap around your wrist and secure it. The wrist strap protects the tonometer from dropping down accidentally. Insert batteries into the tonometer (page 9).

1. Press the measurement button to turn the device on.
2. Following a brief pause, the display will show "LoAd" reminding the user to load the single-use probe into the tonometer prior to measurement.

Loading the probe

1. Accurate measurement is guaranteed only when using probes made by Icare Finland Oy.
2. Open the package, and load the probe into the probe base, according to the illustrations.
3. Place the opened disposable package against the collar component.
4. Turn the device and the probe package at the same time, so that the probe package is uppermost. Check that the probe has dropped into the probe base before removing the package.
5. The probe in the probe base.
6. Raise the device to the vertical operating position ensuring the probe cannot fall out of the probe base. Press the measurement button. The probe will be magnetized by moving back and forward for a short time, and the mechanism that prevents the probe dropping out accidentally will start.

The device is ready for measurement when the reading 00 appears in the display. If the display shows the message LoAd, there is either no probe in the device, or the probe is stuck. Check that the probe has not dropped out, or press the function button again.

Turning the tonometer OFF


1. Press either of the selector buttons until the display shows End.
2. Press the measurement button for two seconds. The display shows ByE and the tonometer turns off. The used probe is ejected from the tonometer. Be careful to dispose of the probe properly (p. 4).
3. The tonometer turns off automatically after two minutes idle time.

Measurement

The measurement should be done without local anesthetic, because a local anesthetic may reduce the tonometer reading*.

Secure the subject and bring the tonometer near to the subject's eye.

The central groove should be in a horizontal position. The distance from the eye to the front part of the collar is the length of the collar. The distance should be 4-8 mm (1/6 - 1/3 inch) from the tip of the probe to the cornea of the eye.

Measure takes place by lightly pressing the measurement button. The tip of the probe should hit the central cornea. Six measurements are made consecutively. Press the measurement button carefully, to avoid shaking the tonometer. After each successful measurement there is a short beep. After the six measurements, there is a longer beep and the IOP is shown on the display after the d (dog/cat) or h (horse) letter.

If there is an erroneous measurement, the tonometer will beep twice and show an error message (page 10). Press the measurement button to clear the error message. If several erroneous measurements appear, see the section "Error messages" (page 10).

Display after the measurements

| Before | After the second measurement | After the sixth measurement |
|--------|------------------------------|-----------------------------|
| 00 | 2.13 | d 13 |

After the sixth measurement, there will be a long beep and the letter d or h appears in the display, followed by the IOP (Intraocular pressure) reading. If the d or h is blinking, it means that the standard deviation of the measurements is greater than normal.

d_ (line down) The standard deviation of different measurements is slightly bigger value than normally, but the effect to the result is not likely to be relevant.

d- (line in the middle) The standard deviation of different measurements is clearly greater than normal, but the effect on the result is usually not relevant. A new measurement is recommended if the IOP is higher than normal IOP.

d° (line up) The standard deviation of different measurements is great and a new measurement is recommended.

*Badouin C, Gastaud P. Influence of topical anesthesia on tonometric values of intraocular pressure. *Ophthalmologica* 1994;208:309-313

Accessing old measurement values

1. Starting position.
2. Press the right or left selector button until Old appears on the display. Then press the measurement button. The old values can now be scanned back and forward by pressing the selector buttons (right= older, left=recent).
3. Press the measurement button to exit the old values search. Old is on the display. Press either of the selector buttons to access the other functions (00=measurement, End=turning OFF).

Changing the calibration table for different species

When the display shows End (p. 4), press first measurement, keep it pressed down and immediately press right selector button. Buttons have to be pressed down at the same time. Then the display shows CAL (or software version). Toggle between these using selector button and choose CAL by pressing the measurement button. After it shows the menu where to choose the calibration table:

do = dog /cat


ho = horse

P = not defined, reserved for future use.

Toggle between these using selector buttons and choose the needed calibration table by pressing the measurement button. When measuring, the result display shows which calibration table is in use (d, H or P). The default selection is d (dog/ cat).

Replacing the probe base


After several months of use, dust may collect in the probe base, affecting the probe movement. The probe base should be replaced if the probe no longer moves smoothly.

1. Unscrew the probe base collar and put it in a safe place.
2. Remove the probe base by tilting the tonometer downward and pull the probe base out of the tonometer.
3. Put a new probe base into the tonometer.
4. Screw the collar in, to lock the probe base.

Replacing the batteries

1. Unscrew the battery compartment locking screw with a screwdriver or a small coin.
2. Remove the battery compartment cover.
3. Remove the old batteries.
4. Put new AA batteries into the tonometer according to figures inside the battery compartment. Terminals should point down on the display side and up on the measurement side.
5. Replace the battery compartment cover and secure it in place by screwing it lightly with the coin or screwdriver. Be careful not to use excessive force when screwing the cover in place.

Error messages:

Press the central button to clear error messages, after which the measurement can be repeated.

The following messages may appear:

- bAtt the batteries are low. Replace the batteries.
- E 01 The probe did not move at all.
If this error message is repeated, turn the tonometer so that the collar faces down for a short time.
If the error message is still repeated take out the probe and replace with a new one. If the error message continues to appear, replace the probe base.
- E02The probe did not touch the eye.
The measurement was made from too far away or the probe base is dirty.
- E03Probe speed was too low.
The measurement was made from too far away or the tonometer was tilted upwards too much.
- E04Probe speed was too high.
The tonometer was tilted downwards. Be sure the groove is in the horizontal position.
- E05The hit was too "soft".
The probe hit the eyelid.
- E06The hit was too "hard".
The probe hit the opening eyelid or calcification in the cornea.
- E07 "Bad bounce".
The probe did not hit the central cornea.

Service procedures

- Replace the batteries when the <bAtt> message appears.
- Change the probe base if the probe does not move smoothly.
- There are no other service procedures that can be carried out by the user. All other servicing and repairs must be carried out by the manufacturer or certified service facilities.
- The device can be cleaned with a damp cloth containing disinfectant.

Technical information

Type: TV01.

Dimensions: 13 – 32 mm (W) * 45 – 80 mm (H) * 230 mm (L).

Weight: 155 g (without batteries), 250 g (4 x AA batteries).

Power supply: 4 x AA batteries.

Measurement range: 1-99 mmHg.

Accuracy of display: ± 1 .

Measurement accuracy: ± 2 (5-30 mmHg), 10% (30-80 mmHg).

Display unit: Millimetre mercury (mmHg).

The serial number is inside the battery compartment cover.

There are no electrical connections from the tonometer to the patient.

The device has B-type electrical shock protection.

Storing/transportation environment:

Temperature +5 to +40 °C

Rel. humidity 10 to 80 % (without condensation)

ISO 13485:2003 certified


CAUTION -Indicates that important operating instructions are included in this User's and Maintenance Manual.

Spare parts and supplies

- Probes in single-use package.
- Probe base replacement kit.

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